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

Contents

	<u>Page</u>
General View	2
Warning and Safety Information	4
Basic pH Meter Quick Reference	5
Installing and Maintaining Electrodes	7
Standardizing for pH Measurement	9
Using Setup	13
Standardizing for Millivolt Measurement	15
Understanding pH Theory	17
Measuring pH	18
Troubleshooting	19
Meter Specifications	21
Accessories	22
CE Marking	23

3. Rear Panel Connectors

Warning and Safety Information

Do not miss out on the benefits of our full guaranty. Complete the guarantee registration card, indicating the date of installation, and return the card to your Sartorius office or dealer.

For safety and operating reasons, only authorized service technicians may open the Basic Meter PB-10 housing. Therefore, only authorized technicians may repair or perform maintenance on this pH meter. Any tampering with the pH meter or negligent or intentional damage to this equipment will void any warranty claims against the manufacturer.

If liquid gets into the pH meter, unplug it from AC power (mains supply) and have an authorized service technician check the pH meter.

If you do not plan to use this pH meter for a relatively long period, please disconnect it from AC power.

For safety reasons, use this equipment only for the application described in this operation manual. Make sure that the buffers used for standardizing have exactly the same values that are stored.


Basic pH Meter Quick Reference


1. Connect power cable to meter power jack and to AC power source.


2. Connect a glass pH electrode to the input and ATC connectors.


3. Press **Mode** until the display indicates the appropriate measurement mode (pH or mV/relative mV).

4. Standardize the meter using up to three buffers by immersing the electrode in a buffer, stirring, then pressing **Standardize** to enter each buffer.

5. The display shows the current reading in pH, mV, or relative mV units.

6. Press **Setup** to review electrode calibration and to clear or select buffer sets.

Installing and Maintaining Electrodes

1. Remove the protective end cover from electrode.

2. Before first use of your glass pH electrode, or whenever the electrode is dry, soak over night in an Electrode Filling Solution, KC1 solution or Electrode Storage Solution.

3. Remove the shorting cap on the BNC connector. Install the combination glass pH/ATC electrode by plugging it into the input connection (push on and twist to lock) and the ATC connector into the ATC jack.

4. Option: Install ORP or Ion Selective Electrode pairs by removing the BNC shorting cap and plugging the BNC connector (twist-lock) into the BNC jack.

5. Rinse and blot-dry electrodes between each measurement (do not wipe). Rinse electrodes with distilled water or deionized water, or part of the next solution to be measured.


6. Store glass pH electrodes in Electrode Filling Solution, KCl solution or Electrode Storage Solution. Always leave the filling hole open during measurement and closed during storage. Refill with Filling Solution when the internal solution level gets low.

Standardizing for pH Measurement


Because electrodes vary in their response, you must standardize your pH meter and electrode to compensate for electrode variation. The more frequently you standardize, the more accurate your measurements. Standardize daily, or more often, for accurate results.

This pH meter allows automatic standardization using up to three buffers.


The pH meter performs automatic temperature compensation.


1. Immerse electrode in a buffer solution. Stir gently. Allow the electrode to reach a stable value.


2. Press and release the **Mode** button until your digital display indicates pH mode. This button toggles between pH and mV modes.


3. Clear existing buffers when doing a new 2 or 3 point standardization. Use the **Setup** and **Enter** buttons to clear existing buffers and to select a new set of buffers. See page 13–14.


4. Press **Standardize**. The meter recognizes the buffer and flashes a buffer icon. When the signal is stable, or when you press **Enter**, the buffer is entered.


5. The meter displays the %slope of the electrode as 100.0%. On entering a second or third buffer, the meter performs a diagnostic check on the electrode and displays the slope.


6. To enter a second buffer, place the electrode in the second buffer solution, stir, allow time for the electrode to stabilize, and press **Standardize** again. The meter recognizes the buffer and displays the first and second buffer icons.


7. Next, the meter performs a diagnostic test of the electrode. The display indicates either *OK* or *Error*. The meter displays the %slope of the electrode.

8. *Error* indicates that your electrode is not working properly. The electrode response must be between 90 and 105% slope. Measurements causing *Error* are not accepted, used or stored by the meter (see "Troubleshooting" on page 19). Press **Enter** to clear the *Error*, then try re-entering the buffer as described in step 6.

9. To set a third standard, place the electrode in the third buffer solution, stir, allow to stabilize, and press **Standardize**. The results will be the same as in steps 7 and 8, except the display will show three buffer icons.

10. After entering each buffer, the meter returns to Measuring operation.

11. Standardize your meter and electrode using at least two buffers with pH values bracketing the expected pH of your samples. Stirring with a magnetic stir bar and stirrer provides faster electrode response.

Note:

If an ATC probe is used, the meter continually adjusts for temperature. Therefore, buffers may vary slightly from the nominal values because of temperature. Default temperature is 25°C.


Note:

The temperature is only displayed when an ATC probe is attached.


Using Setup

The **Setup** button lets you clear all the buffers that you have entered, review calibration information, or select the buffer set that you want.


Note: You can escape the setup mode at any time by pressing **Mode**.


1. Press **Setup** and the meter displays a flashing *Clear* icon. Use this step only when you wish to clear all buffers you have entered. To clear all existing buffers, press **Enter**. The meter clears all buffers and returns to Measuring.


2. Press **Setup** again to show electrode performance. If the meter has accepted two buffers, it will display *OK*, display the slope between the first and second buffers and display the two buffer icons.


3. Pressing **Setup** again shows the electrode slope between the second and third buffers (if three buffers have been entered) and shows the second and third buffer icons.

4. Press **Setup** again to display a flashing *Set* icon and to display the first buffer set icons.

5. Press **Enter** to select the set of buffers shown on the display or Press **Setup** again to view the next set of buffers. Continue pressing **Setup** to view the third buffer sets.

6. Press **Enter** to select the displayed buffer set that contains the buffer you want to use. Press **Setup** again, or press **Mode** at any time, to return to Measuring.


Note: You may mix buffers from different sets.

Standardizing for Millivolt Measurement (Relative Millivolts)


You use millivolt measurement for determining ion concentration and for measuring redox potential (also called ORP, oxidation reduction potential). You will normally use an ion selective electrode (ISE) to measure ion concentration. The ISE senses the ion concentration and responds with a millivolt potential. The millivolt readings are then used to determine ion concentrations.

ORP measurements indicate the oxidizing or reducing capability of a solution. You can use ORP values to monitor or control solutions requiring a set amount of oxidants or reductants.

1. Immerse electrode in a standard solution.


2. Press the **Mode** button until your display indicates *mV* mode.



3. Press **Standardize** to enter a mV standard and read relative mV.

4. When the signal becomes stable, or when you press **Enter**, the current absolute mV value becomes zero relative millivolts.

5. To clear a mV offset and return to absolute millivolt mode, press **Setup**. The meter displays a flashing *Clear* icon, and shows the current relative millivolt offset.

6. To clear the mV standard, press **Enter**. You then return to absolute mV mode.

Understanding pH Theory

Defining pH

The measurement of pH plays an important role in identifying and controlling acidity and alkalinity levels for industry and research. pH is a measure of the acidity or alkalinity of a solution and can be represented by this equation:

$$\text{pH} = -\log [\text{H}^+]$$

with $[\text{H}^+]$ representing the concentration of hydrogen ions in the solution. pH is sometimes referred to as the power of the hydrogen ion in a solution.

By using a pH meter, you can determine exact pH levels of solutions. For example, rather than say that lemon juice is quite acidic, you can say that lemon juice has a pH of 2.4. An exact pH value can be used to control or measure acidity levels for manufacturing processes or for basic research. pH values generally range from 0 to 14, with a pH value of 7 being the neutral point, or the value of pure water. pH values greater than 7 represent increasing alkalinity, whereas pH values below 7 represent increasing acidity (Figure 1).



Figure 1.
pH Scale showing the relative acidity or basicity of some common substances.

Measuring pH


To measure pH with a conventional glass pH electrode, the meter uses a glass bulb that is sensitive to hydrogen ions. The potential developed at the glass membrane is directly related to the pH of the solution.

The glass electrode is paired with a reference electrode, which completes the electrical measuring circuit and provides a stable reference point. These two electrodes are joined to create a combination electrode. The combination glass electrode is connected to the pH meter, which reads the voltage, converts it to pH units, and displays the result.


Combination Glass pH Electrode


Troubleshooting

1. If the signal from the electrode is out of range, the display will show "---". This may happen when the electrode is not in a solution.

2. The meter will display *Error* when it detects an error in electrode response. During standardization, the message indicates that the electrode is less than 90% or more than 105% of the correct response. The *Error* message can indicate either a bad electrode or bad buffer(s).

3. If the meter detects an error in the temperature probe, the display shows ---°C. If no temperature probe is connected the meter uses the default temperature of 25° C.

Electrode Test	
pH 7	0 ± 30 mV
pH 4	159 to 186 mV more than pH 7
pH10	159 to 186 mV less than pH 7

4. To test the pH electrode, place it in a good pH 7 buffer. Press **Mode** to use the mV mode, and note the millivolt reading. Make sure the meter is in mV mode and not relative mV mode. Repeat for either a pH 4 or pH 10 buffer. The electrode signal must be within the limits shown below (when temperature is near 25°C).

5. To test the meter for correct operation, install the BNC (input) shorting cap. Press **Mode** to select the mV mode, and note the mV reading. Make sure the meter is in mV mode and not relative mV mode. If the meter reads 0 ± 0.3 mV, it is measuring correctly. Note that a long-term drift of 0.1 mV / month since last calibration is specified.

Meter Specifications

pH	Range	0 to 14.00
	Resolution	0.01
	Accuracy	± 0.01
mV	Range	± 1500.0 mV
	Resolution	0.1 mV
	Accuracy	± 0.4 mV
Temperature	Range	-5.0 to 105.0°C
	Resolution	0.1°C
	Accuracy	± 0.2 °C
Standardization		0, 1, 2 or 3 buffers
Auto buffer recognition	16 buffers	
		2, 4, 7, 10, 12
		1, 3, 6, 8, 10, 13
		1.68, 4.01, 6.86, 9.18, 12.46
<hr/>		
Auto Temperature Compensation (ATC)		
Automatic electrode slope correction for 90–105%		
Direct reading with a glass electrode		

Accessories

You can order the following accessories for your pH meter:

Plastic-body pH/ATC Electrode, with built-in temperature sensor 3mol/l KC1 liquid-filled	PY-P10
Plastic-body pH/ATC Electrode, with built-in temperature sensor gel-filled	PY-ASI
Glass body pH/ATC Electrode, platinum junction, 3mol/l KC1 liquid-filled	PY-P11
Plastic body pH Electrode, gel-filled	PY-P20
Glass body pH Electrode, platinum junction	PY-P21
ATC Temperature Probe	PY-T01

Other electrodes, such as ion selective electrodes and redox electrodes are also available. Call your sales representative.

CE Marking

The CE marking affixed to the equipment indicates that the equipment meets the requirements of the following Directive(s):

Council Directive 89/336/EEC "Electromagnetic compatibility (EMC)"

Applicable European Standards:

Limitation of emissions:

EN 50081-2 Industrial environment

Defined immunity to interference:

EN 50082-1 Residential, commercial and light industry

EN 50082-2 Industrial environment

Important Note:

The operator shall be responsible for any modifications to Sartorius equipment and for any connections of cables or equipment not supplied by Sartorius and must check and, if necessary, correct these modifications and connections.

On request, Sartorius will provide information on the minimum operating specifications (in accordance with the Standards listed above for defined immunity to interference).

Council Directive 73/23/EEC

"Electrical equipment designed for use within certain voltage limits"

Applicable European Standards:

EN 60950 Safety of information technology equipment including electrical business equipment

EN 61010 Safety requirements for electrical equipment for measurement, control and laboratory use

Part 1: General requirements

If you use electrical equipment in installations and under ambient conditions requiring higher safety standards, you must comply with the provisions as specified in the applicable regulations for installation in your country.